

BS-1: Universal Pointed Tip



Universal pointed tip

The BS-1 is a pointed universal piezo tip that is best for general deposit removal in all quadrants.

BS-2: Flat Chisel Tip



Flat chisel tip

The BS-2 is a flat chisel piezo tip that is best for removing resin-free ortho cements. It may also be used for supragingival calculus removal.

BS-3: Rounded Flat Tip



Rounded flat tip

The BS-3 tip is a rounded flat piezo tip. The rounded edge is good for gentle treatment.

BS-P Piezo Tip



BS-P slim perio tip

The BS-P is a slim perio tip designed for supra and sub-gingival scaling.

BS-PFU Piezo Tip



BS-PFU fine straight perio tip

The BS-PFU is a fine straight-angled perio universal tip. This tip was designed to reach difficult-to-reach spots near the posterior teeth.

BS-PFL Piezo Tip



BS-PFL fine left-angled perio tip

The BS-PFL tip is a fine left-angled tip designed for perio procedures and for reaching difficult-to-reach spots near the posterior teeth.

BS-PFR Piezo Tip



BS-PFR Fine right-angled perio tip

The BS-PFR is a fine right-angled perio tip. This tip was designed to reach difficult-to-reach spots near the posterior teeth.

BS-10Z Piezo Tip



BS-10Z Wide flat-headed perio tip

The BS-10Z is a wide, flat headed, straight-angled tip specially designed for perio functions and for accessing difficult-to-reach spots near the posterior teeth. This tip is recommended for work on medium pockets less than 4mm deep.

BS-H2R Piezo Tip



BS-H2R Curved right-angled tip

The BS-H2R is a curved right-angled tip specially designed for perio functions and for accessing difficult-to-reach spots near the posterior teeth. This tip is recommended for treatment on premolars, molars, furcation, and narrow spaces.

BS-H2L Piezo Tip



BS-H2L Curved left-angled tip

The BS-H2L is a curved left-angled tip specially designed for perio functions and for reaching difficult-to-reach spots near the posterior teeth. This tip is recommended for treatment on premolars, molars, furcation, and narrow spaces.

BS-10P Piezo Tip



BS-10P 45° angled tip

The BS-10P is a 45° angled tip recommended for procedures on shallow pockets.